

INSTRUCTIONS FOR USE

Ehrlichia chaffeensis and Anaplasma phagocytophilum IFA IgM Antibody Kit

Catalog Number: E21M-120
Size: 120 test
Storage: 2 – 8 °C

An Indirect immunofluorescence assay for the detection and semi-quantitative determination of IgM class antibody against both *Ehrlichia chaffeensis* and *Anaplasma phagocytophilum* in human serum or plasma

For in-vitro diagnostic use only



1135 E. Truslow Avenue
Fullerton, California 92831 USA
Phone: +1-714-525-7660
Fax +1-714-525-7614
Email: info@fullerlabs.com
URL: www.fullerlabs.com



MediMark Europe Sarl
11, rue Émile Zola – BP 2332
F-38033 Grenoble Cedex 2 – France

INTENDED USE

The *Ehrlichia chaffeensis* and *Anaplasma phagocytophilum* IFA Human IgM Antibody kit is intended for the detection and semi-quantitation of IgM class human antibody to both ***Ehrlichia chaffeensis*** and ***Anaplasma phagocytophilum*** by micro-immunofluorescence assay.

SUMMARY AND EXPLANATION OF TEST

Both *Ehrlichia chaffeensis* and *Anaplasma phagocytophilum* are tick-borne human pathogens. The microimmunofluorescence assay utilizes semi-purified elementary bodies and morulae from cell culture-propagated organisms. For optimal adhesion and background contrast, these antigens are dispersed in a matrix of sonicated cellular material.


Patient sera are pre-treated to remove sources of rheumatoid factor, then diluted in buffered saline. They are then incubated in the individual slide wells to allow reaction of patient antibody with the solid-phase antigens. Slides are then washed to remove unreacted serum proteins, and fluorescence-labeled anti-human IgM (Conjugate) is added. This Conjugate is allowed time to react with antigen-antibody complexes. Then removed via another wash step. The resulting reactions can be visualized using standard fluorescence microscopy, where a positive reaction is seen as sharply defined apple-green fluorescent inclusions (morulae) and small cocci within a red cellular matrix background. As the slides are viewed through the microscope from left-to-right (frosted-end to the left), the *A. phagocytophilum* antigen is viewed first. To the right of this antigen, in each slide well, is the *E. chaffeensis* antigen dot. Positive reactions at screening dilution should be retested at higher dilutions to determine the highest reactive or endpoint dilution.

REAGENTS AND MATERIALS SUPPLIED

IFA Ag x 12

Substrate Slides (10)

10 x 12-well masked slides containing antigen dots of (LEFT) *Anaplasma phagocytophilum* and (RIGHT) *Ehrlichia chaffeensis*, both antigens dispersed a cellular background matrix. Slides are fixed and packaged under vacuum, ready to use.



CONJ FITC

IgM Conjugate, 2.5 mL

Dropper bottle with a yellow cap contains affinity-purified donkey anti-human IgM (heavy chain-specific) conjugated to the fluorophore DyLight 488, with bovine serum albumin and Evans' blue counterstain.

SAMP DIL

IgM Sample Diluent, 15 mL

Protein-based buffer containing goat anti-human IgG, ready for use.

CONT +

EC Positive Control, 0.5 mL

Dropper bottle with a blue cap contains *E. chaffeensis*-reactive serum pre-treated and provided at a 1:64 screening dilution. Endpoint titer is 1:512 (1:256-1:1024) to the *Ehrlichia* antigen spot.

CONT +**AP Positive Control, 0.5 mL**

Dropper bottle with a blue cap contains A. phagocytophilum-reactive serum pre-treated and provided at a 1:64 screening dilution. Endpoint titer is 1:512 (1:256-1:1024) to the Anaplasma antigen spot.

CONT -**Negative Control, 0.5 mL**

Dropper bottle with a red cap contains non-reactive serum, provided at a 1:64 screening dilution.

MM**Mounting Medium, 1 mL**

Dropper bottle with a white cap contains 50% glycerol in PBS.

BUF WASH PBS**PBS, 1 liter**

Add supplied powder to 1 liter purified water to produce phosphate-buffered saline at pH 7.2. Mix well.

Warnings

1. Since no testing can assure the absence of infectious agents, however, the Controls, as well as all serum specimens and equipment coming in contact with these specimens, should be handled with good laboratory practices to avoid skin contact and ingestion.
2. The substrate slides are prepared with chemically inactivated antigens. However, the slides should be considered potentially infectious and handled accordingly.

Storage and Handling

Kit components should be stored at 2-8°C. Bring them to room temperature (20°-25°C) before opening bottles or slide envelopes.

SPECIMEN COLLECTION

Allow blood samples to clot and separate sera by centrifugation. Transfer sera aseptically to tightly closing sterile containers. Store at 2-8°C. If testing is to be delayed longer than 5 days, freezing the sample at -20°C or colder is recommended. Acute specimens should be drawn at the onset of illness; convalescent specimens should be obtained at two and four week intervals to check for titer changes.

PROCEDURE

The kit supplies sufficient materials for 120 determinations.

Materials Required But Not Supplied

- Distilled or deionized water
- Clean 250 or 500 mL wash bottle for PBS
- Test tubes or microtiter plate for serum dilutions
- Precision pipette(s)
- 24 x 50 mm glass cover slips
- Fluorescence microscope with filter system for FITC (maximum excitation wavelength 490 nm, mean emission wavelength 530 nm) and 400X magnification
- 37° water bath or incubator
- Humid chamber for slide incubation steps

Precautions

- Do not use components past expiration date.
- Conjugate is photosensitive and is packaged in opaque plastic for protection. Store in the dark and return to storage after use.
- Conjugate contains Evans' blue dye, which may be carcinogenic. Avoid contact with skin.
- Liquid reagents contain thimerosal at 0.01%, which may be toxic if ingested.

PREPARATION OF SAMPLES AND REAGENTS

1. **Prepare Wash Buffer** by adding contents of PBS packet to 1 liter purified water and mixing thoroughly:
2. **Prepare screening dilutions** of patient sera by making an initial 1:16 dilution using Sample Diluent in micro-centrifuge tubes. Mix and allow a minimum of 5 minutes for the reaction, then centrifuge at high speed to remove the aggregated IgG. Dilute 10 µL of this supernate with 30 µL Wash Buffer, resulting in a final 1:64 screening dilution.

ASSAY PROCEDURE

Allow all reagents and sera to reach ambient temperature before starting timed assay procedure.

1. Prepare dilutions of each Positive Control to include 1 dilution above the stated endpoint and one dilution below (i.e. 1:256-1:1024). These Controls are bottled at 1:64.
2. For each serum screening dilution prepared (see above), add 10 µL to a slide well and record the location for later reference. For each assay run, include the Negative Control and dilutions of the Positive Controls prepared above.
3. Place slides in a humid chamber and incubate for 30 minutes at 37°± 0.5°C.
4. Remove humid chamber from incubator. Also remove conjugate from storage. Rinse slide wells with gentle stream of PBS from wash bottle. Shake or tap beaded PBS from slides into a sink, then repeat this wash step 3X without allowing the wells to dry.
5. To each slide well, add 1 drop (10-15 µL) Conjugate and then return slides to the humid chamber for 30 minutes incubation at 37°± 0.5°C. Incubation should be in the dark to protect the photosensitive Conjugate.
6. Wash slides as in step 5, above.
7. Add 2-3 drops Mounting Medium to each slide and carefully apply cover glass.
8. Read the stained substrate slides at 400X magnification, comparing each well to the visual intensity, antigen density and appearance of Positive and Negative Control wells. Slides may be stored at 2-8°C in the dark for up to 24 hours.

QUALITY CONTROL

The Negative Control and dilutions of the Positive Controls should be assayed with each daily run. The Negative Control well is an example of a non-reactive serum, with either uniform red counterstain or slight, but uniform greenish staining. The Positive Control wells should give endpoint titers from 1:256 to 1:1024 with the appropriate Positive Control. The fluorescence intensity at 1:512 may be used as the cut-off level required for a patient reaction to be called positive. If either of the Controls does not react as specified, the assay run should be considered void, reagent components and procedural steps should be rechecked, and the assay repeated from step #1.

The Negative Control well is an example of fluorescence patterns that are to be considered negative. If characteristic morulae and elementary bodies are seen in this well, similar to that seen in the Positive Control wells, there has been a breakdown in technique and the assay must be repeated.

INTERPRETATION OF RESULTS

A positive reaction appears as fluorescent particles within a red background matrix. The size, appearance and density of the reaction must be compared with the Positive and Negative Control reactions.

Patient Specimens

Negative at 1:64: Report as negative for the respective antibody specificity. Further serum specimens should be drawn if the original was taken soon after onset and ehrlichiosis is still suspected.

Positive at 1:64 and greater: Serum titers at 1:80-1:320 suggest 1) titers preceding peak levels (early), 2) titers after peak levels (past exposure) or 3) titers reflecting cross-reactivity to a related organism (ie. *Ehrlichia* spp). Titers greater than 1:320 and/or IgM titers, when present, are a reliable indicator of recent infection.

Paired Sera: A four-fold increase in titer between acute and convalescent sera supports the diagnosis of recent infection by *Ehrlichia chaffeensis*, *Anaplasma phagocytophilum* or a closely related organism.

LIMITATIONS

Crossreaction between *Ehrlichia chaffeensis* and both *Ehrlichia canis* and *Ehrlichia ewingii* by IFA is variable from moderate to strong, and can be differentiated by a variety of alternate methods including western immunoblot technique.

SPECIFIC PERFORMANCE CHARACTERISTICS

Anaplasma specificity was tested by 95 sera from a non-endemic region. All 95 of these sera had titers <1:80. Twelve sera from a regional public reference laboratory were also tested for concordance. All 8 positive sera were detected with titers within 1 dilution and the 4 negative sera were each <1:80. For sera with *E. chaffeensis* IFA titers >1:80 by standard IFA protocols were all positive in this assay, as well.

Sera from *E. chaffeensis* non-endemic regions were tested in-house, 159 from New York and 120 from Southern California. There were no positives. (100% specific). As the New York sera were from an endemic region for *Anaplasma phagocytophilum*, there were found 14 dogs (8.8%) seropositive for this related organism.

REFERENCES

1. Dumler, J. S., K. M. Asanovich, J. S. Bakken, P. Richter, R. Kimsey, and J. E. Madigan. 1995. Serologic cross-reactions among *Ehrlichia equi*, *Ehrlichia phagocytophila*, and human granulocytic *Ehrlichia*. J. Clin. Microbiol. 33:1098-1103.
2. IJdo, J. W., Y. Zhang, E. Hodzic, L. A. Magnarelli, M. L. Wilson, S. R. Telford III, S. W. Barthold, and E. Fikrig. 1997. The early humoral response in human granulocytic ehrlichiosis. J. Infect. Dis. 176:687-692.
3. Magnarelli, L. A., J. W. IJdo, J. S. Dumler, R. Heimer, and E. Fikrig. 1998. Reactivity of human sera to different strains of granulocytic ehrlichiosis in immunodiagnostic assays. J. Infect. Dis. 178:1835-1838.
4. Walls, J. L., M. Aguero-Rosenfeld, J. S. Bakken, J. L. Goodman, D. Hossain, R. C. Johnson, and J. S. Dumler. 1999. Inter- and intralaboratory comparison of *Ehrlichia equi* and human granulocytic ehrlichiosis (HGE) agent strains for serodiagnosis of HGE by the immunofluorescent-antibody test. J. Clin. Microbiol. 37:2968-2973.
5. Ristic, M., D.L. Huxsoll, Weisiger, R.M., Hildebrandt, P.K., and Nyindo, M.B.A. 1972. Serological Diagnosis of tropical Canine Pancytopenia by Indirect Immunofluorescence. Infect. Immun. 6:226-231.
6. McBride, J.W., Corstvet, R.E., Gaunt, S.D., Boudreaux, C, Guedry, T, and Walker, D.H. 2003. Kinetics of Antibody Response to *Ehrlichia canis* Immunoreactive proteins. Infect. Immun. 71:2516-2524.
7. Ohashi, N., Unver, A., Zhi, N., and Rikihisa, Y. 1998. Cloning and Characterization of Multigenes Encoding the Immunodominant 30-Kilodalton Major outer membrane Proteins of *Ehrlichia canis* and Application to the Recombinant Protein for Serodiagnosis. J. Clin. Microbiol. 36:2671-2680.
8. Unver, A., Rikihisa, Y., Ohashi, N., Cullman, L., Buller, R., and Storch, G.A. Western and Dot Blotting Analyses of *Ehrlichia chaffeensis* Indirect Fluorescent-Antibody Assay-Positive and -Negative Human Sera by Using native and Recombinant *E. chaffeensis* and *E. canis* Antigens. J. Clin. Microbiol. 37:3888-3895.
9. Breitschwerdt, E.B., Hegarty, B.C. and Hancock, S.I. 1998. Sequential Evaluation of Dogs Naturally Infected with *Ehrlichia canis*, *Ehrlichia chaffeensis*, *Ehrlichia ewingii*, or *Bartonella vinsonii*. J. Clin. Microbiol. 36:2645-2651.

Original 7/95

Current Version D (4/05)