# Immunodiagnostic Assays for Flea-Borne Typhus

Elisabeth Laderman, Oscar Aurelio, Julia Reyna and Lee Fuller

## Abstract

Flea-borne typhus in humans is found worldwide and is etiologically related to both *Rickettsia typhi* and *Rickettsia felis*. Endemic foci are found in many regions, including Los Angeles and Orange Counties in Southern California. This poster describes some of our testing with both a 2-antigen micro-IFA and specific ELISA assays for these two agents. Although multiple assays for *R. typhi* have been commercially available for over 20 years, the *R. felis* ELISA has been developed more recently. Native rOmpB from octyl-glucoside-extracted *R. felis* was eluted in a pure form in the void volume of Sephadex G-100 gel chromatography. This assay is compared with values derived using the 2-antigen micro-IFA assay using normal healthy blood donor sera from the Texas Gulf Coast.

## Materials and Methods

**Purification of rOmpB from *R. felis*:** *R. felis* was propagated in D.mel-2 cells adapted to serum-free medium at 28°C. Harvests were pooled and stored at -85°C. Rickettsia were harvested by differential centrifugation, purified via standard Renografin density gradient and stored at -85°C. An aliquot of *R. felis* in PBS was then brought to 1% octyl- $\beta$ -D-glucopyranoside and extracted for 15 minutes at 25°C. The cell-free supernate was removed after centrifugation and layered over Sephadex G-100 equilibrated with PBS. Fractions were collected by gravity flow and analyzed by a BCA protein assay (Pierce) and SDS-PAGE using 10% Bis-Tris gels (Life Technologies). A single peak of 120 kDa protein was found in the void volume and these peak fractions were pooled and stored in aliquots at -85°C.

**MIF assay for flea-borne typhus:** *R. felis* was propagated as described (above) and *R. typhi* was propagated in Vero-6 cells at 32°C. Both were purified as described and inactivated for slide preparation in 0.25% methanol-free formaldehyde. Antigens were blended to a standard density in a background matrix, robotically applied to 12-well teflon-masked IFA slide stock, acetone-fixed and packaged under vacuum.

***R. felis* ELISA development:** Purified rOmpB isolated as described above was coated onto medium-binding ELISA wells (Costar) in a series of dilutions. Sera which were negative, weakly or strongly positive for *R. felis* by MIF assay (above) were used to optimize the titer of antigen and horseradish peroxidase (HRP) conjugated anti-human IgG. Test sera were diluted 1:100 in sample diluent and added to the coated wells. A serum incubation of 60 minutes and conjugate incubation of 30 minutes (both at ambient temperature) were separated by 4X washing with PBST. TMB substrate was added for a 10 minute incubation followed by a stop solution. Results were measured at 450 nm on an ELISA plate reader and compared with MIF titers, especially with sera giving weak reactions at 1:32 (typhi) or 1:64 (felis) using the MIF technique. These sera define the cut-off calibrator for the ELISA results, with absorbance values within 10% of this cutoff absorbance value defining an “equivocal” range. Titers within this equivocal range are expected to be evenly divided between high negative and low positive.

***R. typhi* ELISA formats:** Two formats were used from Inventory, with kit components including ELISA modules coated with either rOmpB antigen or LPS antigen.

# Individual Serum Data

180 serum samples (n= 180) from the Gulf Coast Regional Blood Center were randomly selected from plasma units processed May 13-14, 2013 that had been cleared by required donor blood testing. Donor age ranged from 17-75 years of age and donor residence is expected to be within the region serviced by Gulf Coast Regional Blood Center in Houston, Texas.

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**Table 1.** Seroprevalence of *R. felis* and *R. typhi* in a population of normal healthy blood donors in an endemic region. The Equivocal range was initially defined as MIF titers between 1:32-1:64 for each ELISA assay, with the ELISA Cutoff Calibrator defined by these minimal positive sera. For total RT prevalence values positive reactivity with either rOmpB or LPS antigens was counted as a positive serum. The correlation of these assays demonstrates the higher sensitivity with the LPS antigen and the LPS-centric sensitivity of the MIF assay. For RF assays the lack of competing antigens allows for better correlation.

TABLE 2

<i>Rickettsia felis</i>					
Result		POS	EQ	NEG	Total
<i>R. typhi</i>	POS	3	2	11	16
	EQ	1	1	11	13
	NEG	13	10	128	151
	Total	17	13	150	180

**Table 2. Comparison of *R. felis* and *R. typhi* results in normal healthy blood donors.** Results were classified as Positive if at least 1 assay result was positive. The Equivocal designation is for sera with at least 1 equivocal titer and no positive results. Only 3 sera were positive for both agents (1.7%), with another 4 sera (2.2%) showing equivocal reactivity to both. The results with the 2 rOmpB ELISA assays show no dual reactivity (0%), further suggesting that dual-reactivity seen in endemic typhus is more probably a true dual reactivity based on multiple flea bites.

## Conclusions

- For serodiagnosis of *R. felis* infection the MIF and ELISA assays correlate closely, even with the low reactivity of healthy blood donor sera.
- For serodiagnosis of *R. typhi* infection both rOmpB and lipopolysaccharide (LPS) antigens are important, with 50% of positive sera reacting with only a single antigen.
- Our data suggests that the *R. typhi* MIF antigen we used is much more sensitive to anti-LPS reactivity.
- Seropositivity rates are similar for *R. felis* and *R. typhi*.
- Our data suggests that these assays could all be screened at lower dilutions or lower cutoffs without losing specificity.